



HB2151 菌株

菌株名称: E. coli HB2151

平台编号: bio-108787

规格: 2ml 甘油管

状态: 20%甘油保存

性质参数: K12 (lac-pro), ara, nalr, thi/F'[proAB, lacIq, lacZ M15].

建议:

1. 菌株保存在 20%甘油的培养基中, 接到菌株后, 请挑取菌液在 2x YT-N 平板 (2xYT medium containing 100 µg/ml nalidixic acid) 上划线, 37°C过夜培养后,

再次挑取单个菌落在 Minimal Medium 平板上划线, 37°C过夜培养后即可使用, 以便使纤毛充分表达, 利于噬菌体感染;

2. 在 Minimal Medium 平板上生长的细菌, 不宜在冰箱中放置过久再使用;

3. Minimal Medium plate 配制方法: Prepare stocks of the following:

1 M MgCl₂•6H₂O: Dissolve 20.33 g in distilled water to a final volume of 100 ml and autoclave.

1 M CaCl₂•2H₂O: Dissolve 14.7 g in distilled water to a final volume of 100 ml and autoclave.

1 M thiamine hydrochloride: Dissolve 33.73 g in distilled water to a final volume of 100 ml and sterilize using a 0.22 µm filter.

20% glucose: Dissolve 20 g of D-(+)-glucose (anhydrous) in distilled water to a final volume of 100 ml and sterilize by filtration through a 0.2 µm filter. Do not autoclave.

In a 500 ml bottle, dissolve 6 g of Na₂HPO₄ (dibasic), 3 g of KH₂PO₄(monobasic) and 1 g of NH₄Cl in distilled water to a final volume of 500 ml. adjust the pH to 7.4 with NaOH.

In a separate 1 liter bottle, add distilled water to 15 g of Bacto-agar to a final Volume of 500 ml. Autoclave both bottles simultaneously to sterilize. Cool both bottles to 50-60°C and combine. Add 1 ml of 1 M MgCl₂•6H₂O, 1 ml of 1 M CaCl₂•2H₂O, 1 ml of 1 M thiamine hydrochloride and 5 ml of 20% glucose. Pour plates immediately.